

Belscan Continental BV Stokkelaar 4

9160 Lokeren

Belgium

Tel: +32 9 340 56 26

Email: info@belscan.be Home: www.Belscan.be

PRODUCT DATA SHEET						
	XS	S	М	L	XL	XXL
Our reference	9002115	9001263	9001264	9001265	9001266	9001267
Description	Nitrile gloves		HY@PRO	N M4	AQL	1,5
	Nitrile gloves black powder-free M4 100 pcs/dispenser box HY@PRO					

PRODUCT PACKAGING IMAGE



PRODUCT DESCRIPTION										
Material			Nitrile							
Colour			Black							
Characteristics		Powderfree		Non Sterile	Ambidextrous					
Surface		Finger Textured								

SIZES										
	XS (5-6)	S (6-7)	M (7-8)	L (8-9)	XL (9-10)	XXL (10-11)	XXXL (11-12)			
Width	70 ± 10 mm	80 ± 10 mm	95 ± 10 mm	110 ± 10 mm	120 ± 10 mm	130 ± 10 mm				
length	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm				

THICKNESS										
	Cuff	Palm	Fingertips	Overal weight (M)						
Thickness	0,05 mm	0,07 mm	0,08 mm	4,0 gram						

PRODUCT DATA SHEET						
	XS	S	М	L	XL	XXL
Our reference	9002115	9001263	9001264	9001265	9001266	9001267
Description	Nitrile gloves					
	Nitrile gloves 100 pcs/dispe	· ·		4		Y


STANDARDI	SATION EN388	3						
EN 388 Mechanical risks	Abrasion resistance	Blade cut resistance	Tear resistance	Punture resistance	Blad cut resistance TDM-test	Impact test		EN 388:2016
Level								
STANDARDI	SATION EN374							
EN 374-1	The gloves must not	change their shape to	oo much, nor must th	ey break or stick.			~	
EN 374-2	The material must no	ot contain any perfora	ations				~	
EN 374-3	Chemicals must not	nemicals must not pass through the glove material						
EN 374-4	Tests are carried out	Tests are carried out on the duration and the amount of chemical substances that penetrate the glove.						
EN 374-4	Chemical		Code letter	Permeation level		Degradation		
	Sodium Hydroxide		K	6		-38,40%	~	
	Hydrogen Peroxide		Р	2		17,60%	~	
	n-heptane		J			27,40%	~	
	Ammonium hydroxid	de	0			0,299	~	
	Formaldehyde		Т	5		0,466	~	
Accelerators							\circ	
EN 374-5 Microorganism	The glove is tight aga	ainst microorganisms	(viral, bacteria and fu	ngi).				EN ISO 374-8 2216
tightness	Testing method		•	VIRUS				
STANDARDI	SATION EN165	523						
EN16523-1		aterial resistance to pe ditions of continuous		als - Part 1: Permeatio	on by potentially hazai	rdous liquid	~	

STANDARDIS	STANDARDISATION EN455									
EN 455 medical gloves for single use	The glove meets the requirements according to EN 455-1, EN 455-2, EN 455-3, EN 455-4	~	EN455							
EN 455-1	The glove has an AQL $<$ 1.5 in regards to the water retention test (sampling inspection in acc. to ISO 2859-1, general Inspection Level I)	~	EN455-1							
EN 455-2:	The tensile strength, thickness and fit is tested with a minimum standard of 6 Newtons.									
EN 455-3	This biological evaluation tests for the presence of harmful substances such as powder and micro-organisms.	~	EN455-3							
EN 455-4	The product tested shelf life	~	EN455-4							
STANDARDISATION EN21420										
ISO EN 21420 Protective gloves	The glove meets the requirements according to EN 420 / ISO EN 21420	~								

Rev-nr: 001 Date: 27/06/2024

PRODUCT DATA SHEET						
	XS	S	М	L	XL	XXL
Our reference	9002115	9001263	9001264	9001265	9001266	9001267
Description	Nitrile gloves					
	Nitrile gloves 100 pcs/dispe	•		4		Y

REGALUTOR	Y AFFAIRS								
PPE-Regulation (EU) 2016/425			Category III						
MD-Regulation (EU) 2017/745									
CE Regulation		SATRA 2777							
EUROPEAN REPRESENTATIVE							\circ	EC REP	
Food Contact (EG) 1935/2004	acidic foods	aqueous foods	fatty foots	alcoholic foods	dry foods		~	Ω"	
LATEX	Does not contain late	ex					~	DALEA	
Storage conditions	Store in original packaging in a dry and dark place at 10°C to 30°C. Refer to guidelines of storage of rubber product as decribed in ISO2230:2002. Ensure that storage is kept cool, dry and dust-free, avoid ventilation and storage close to phocopy equipment. Copper-ions discolour the glove. Shelf life (months) 36 Protect gloves against ultraviolet light sources, as sunlight and oxidizing agents. Storage above 30°C will lead to accelerated aging and should be avoided.							***	
Single/double/Reus able	Examination glove fo	or single use					~	2	
Sterile / Non Sterile	Non Sterile glove						~	NOR	
ASTM F739	Symbol used for SW	chemotherany tested	I gloves that have bee	n tested either by AS	TM E739 or by ASTM	D6078		TESTED FOR USE WITH	
ASTM D6978	Symbol used for SW	chemotherapy tested	i gloves that have bee	in tested either by As	TIVIT 739 OF BY ASTIVIT	50378		CHEMOTHERAPY DRUGS	
REACH Compliant	REACH is a regulation of the European Union, adopted to improve the protection of human health and the environment from the risks that can be posed by chemicals, while enhancing the competitiveness of the EU chemicals industry.						~		
CE	CE conformity							(€	
Intrastat number	4015190000								
Made in	CN								



LOGISTIC DATA SALES UNIT

General information		Material	Carton					
Size	pieces	EAN Code	REF code	Length (mm)	Width (mm)	Heigth (mm)		Gross weight
XS	100	5404017208131	9002115	22	12	6	350	382
S	100	5412276251935	9001263	22	12	6	385	417
М	100	5412276251959	9001264	22	12	6	400	432
L	100	5412276251973	9001265	22	12	6	460	492
XL	100	5412276251997	9001266	22	12	6	520	552
XXL	80	5412276252017	9001267	22	12	6	560	592

LOGISTIC DATA OUTER PACKAGING

General information		Material	Carton					
Size	pieces	EAN Code	REF code	Length (mm)	Width (mm)	Heigth (mm)		Gross weight
XS	10	5404017208148	9002115	315	250	230	3,82	4,17
S	10	5404017205703	9001263	315	250	230	4,17	4,52
М	10	5404017205727	9001264	315	250	230	4,32	4,67
L	10	5404017205741	9001265	315	250	230	4,92	5,27
XL	10	5404017205765	9001266	315	250	230	5,52	5,87
XXL	10	5412276252024	9001267	315	250	230	5,92	6,27

LOGISTIC DATA PALLET

General information		Pallet: EURO-PALLET							
Size	Cartons per layer	Layers per pallet	Netto height	Gross Heigth (incl. Euro pallet 14,4cm)	Gross weight (incl. Euro pallet 25kg)				
XS	11	6	25	164,4	300,22kg				
S	11	6	25	164,4	323,32kg				
М	11	6	25	164,4	333,22kg				
L	11	6	25	164,4	372,82kg				
XL	11	6	25	164,4	412,42kg				
XXL	11	6	25	164,4	438,82kg				

PRODUCT DATA SHEET						
	XS	S	М	L	XL	XXL
Our reference	9002115	9001263	9001264	9001265	9001266	9001267
Description	Nitrile gloves					
	Nitrile gloves black powder-free M4 100 pcs/dispenser box HY@PRO			4		Y

Declaration of conformity Food



Food Compliance Declaration - Belscan

Date: 04-07-2022 Owner: Q-Team Code: I.02.09.03.DOC.10 Version: 001 Authorizer: QA

Declaration of compliance for materials and articles intended to come into contact with food

The Supplier, or his authorized representative who is established in the community:

Name: Belscan Continental BV

Address: Stokkelaar 4-9160 Lokeren, Belgium

Declares that the following product(s):

Ref. 9001263-67-9002115: Nitrile gloves black powder-free M4 100 pcs/dispenser box HY@PRO

is in accordance with Regulations (EC) No 2023/2006 (Good Manufacturing Practices)¹, and No. 1935/2004 (Non-Plastic Materials)².

For use in the following circumstances

- type of food stuff intended for coming into contact with materials / objects: All types of food
- Storage time and temperature of the material/object: 3 years- between 5° c and 30 °c
- Possible treatment for the material/object: not suitable for heat treatment
- Ratio surface/volume: 10 mg/d²

(Appropriate information on all substances for which there are restrictions, at the level of the EU as well as at the Belgian level, so that all future users can comply with those restrictions. In the absence of any national or European regulation, all information on international restrictions, norms or guide values has to be provided (Council of Europe, WHO, Codex Alimentarius...))

Completed at,

Lokeren

Radea Ahaddour Quality department For Belscan Continental BV



Rev-nr: 001 Date: 27/06/2024